
AeroEclipse II

Breath Actuated Nebulizer (BAN)

Instructions for Use

- Ensure mouthpiece is inserted into side opening of the nebulizer with the exhalation valve facing down. Inspect the tubing for fray, wear or loose particulate, and ensure nothing is blocking the air supply pathway of the nebulizer.
- Prior to use, inhale and exhale through the nebulizer mouthpiece to verify the movement of the green feedback button located on top of the device.
- Unscrew and remove top of nebulizer. Place prescribed medication into nebulizer cup (maximum cup volume 6 mL). Reattach top and gently hand tighten.
- Attach one end of the supplied tubing to the fitting in the bottom of the nebulizer and the other end to the flowmeter or compressor. Ensure both ends are securely engaged.
- Set the flowmeter to 7 to 8 litres per minute (Lpm) with an air source capable of delivering 50 PSI (344.7 kPa).
- Place mouthpiece in mouth, and inhale slowly and deeply. As you inhale, the green feedback button on top of the nebulizer will move into the fully down position, indicating that the **AEROECLIPSE**® II BAN is producing aerosol in response to inhalation. Exhale normally. A valve on the mouthpiece opens allowing passive exhalation, as the green feedback button returns to the up position indicating no aerosol is being produced. Do not place the lips over the exhalation valve on the bottom of the mouthpiece. Doing so will prevent the valve from functioning properly.
- Continue to inhale from the nebulizer until nebulizer begins to sputter, or as instructed by your physician or qualified healthcare worker.

The **AEROECLIPSE**® II BAN can be operated in Continuous Mode by rotating the mode selector 90 degrees clockwise to the continuous mode position (▲). The nebulizer can be returned to the Breath Actuated Mode by rotating the mode selector 90 degrees counterclockwise to the breath actuated position (◡).

Cleaning Instructions

The mode switch should be rotated to the breath actuated position (green button is up, mode switch has rotated counterclockwise). Remove the supply tubing from the bottom of the nebulizer and unscrew and remove the nebulizer top.

Between treatments with the same patient, the nebulizer may be processed according to any of the methods listed below.

- Sterile Water:** Thoroughly rinse the disassembled parts in sterile water and shake to remove excess water.
- Isopropyl Alcohol:** Wipe the outside surfaces with a tissue dipped in 70% Isopropyl Alcohol then submerge all parts in a soak of 70% isopropyl alcohol for 1 minute before shaking to remove any excess alcohol.
- Air Dry:** Disassemble the device and shake out any excess drug.
- Hand Washing:** Wash nebulizer parts in a warm water/liquid dish detergent solution and rinse in clean warm water.

For all methods, reconnect the bottom of the nebulizer to the compressed medical gas supply to clear the orifice. Always allow the parts to air dry thoroughly before reassembling.

Note: Always follow accepted procedures for cleaning between uses with the same patient.

Single Patient Use


The device should be replaced after approximately 7 days or immediately if components are cracked or the nebulizer does not function properly. This device is not intended to be reprocessed and should not be autoclaved or boiled.

Disposal

The **AEROECLIPSE**® II BAN can be disposed of with domestic waste unless this is prohibited by the disposal regulations prevailing in the respective member countries.

CEN 13544-1 Data

The Aerosol Output Rate for the **AEROECLIPSE**® II BAN is 0.17 ± 0.02 mL/min and the Aerosol Output is 0.61 ± 0.04 mL in Breath Actuated Mode.

Indications for Use

The **AEROECLIPSE**® II BAN is intended to be used by patients who are under the care or treatment of a licensed healthcare provider. The device is intended to be used by these patients to administer aerosolized medication prescribed by a physician or healthcare provider. The intended environments for use include the home, hospitals and clinics.

⚠ WARNING

Close supervision is necessary when this product is used by or near infants, children or anyone requiring assistance to use the device. It contains parts that may present a choking or strangulation hazard. DO NOT LEAVE CHILDREN UNATTENDED DURING TREATMENTS.

⚠ CAUTIONS

- The **AEROECLIPSE**® II Breath Actuated Nebulizer can deliver aerosolized medication in either a Breath Actuated or Continuous Mode. In Breath Actuated Mode, medication is delivered on inspiration only and may result in an increase of the dosage to the lungs. Healthcare providers should take this information into consideration with regards to dosage.
- Failure to properly follow the above cleaning instructions may have an adverse affect on device performance.

REORDER NUMBER

AEROECLIPSE® II BAN 10550191050

TRUPELL MEDICAL INTERNATIONAL

725 Baransway Drive,
London, Ontario, Canada N5V 5G4
☎ customerservice@truellmed.com

☎ +1-519-455-7060
☎ +1-519-455-7858

www.truellmed.com

Canadian Hospital Distribution:

Truell Healthcare Solutions
758 Baransway Drive, London, Ontario, Canada N5V 5J7
☎ customerservice@truellhls.com
☎ +1-800-265-5494
www.truellhls.com

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