**DETENDEUR
 NOVA 40
 PRESSURE-REGULATOR**

SOMMAIRE

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CONTENTS

1. ABBREVIATIONS AND SYMBOLS
2. AREA OF APPLICATION
3. FEATURES
4. SAFETY
5. INSTALLATION
6. START-UP
7. MAINTENANCE

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In order to maintain the quality of our product in the best safety conditions throughout its period of use, we strongly recommend you to read this brochure and to follow the instructions contained within.

Ind.	Modifications	Nom	Date	Ind.	Modifications	Name	Date
+	Edition originale	ERO	04/07	+	Original issue	ERO	04/07

52, rue de Lagny
 93100 MONTREUIL - FRANCE
 Tél : +33 (0)1.41.72.90.00
 Fax : +33 (0)1.41.72.90.01
 E-mail : cahouet@cahouet.com
 www.cahouet.com

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
1. ABBREVIATIONS AND SYMBOLS

1/1 Abbreviations

DET = Pressure regulator
P1 = Nominal inlet pressure
P2 = Nominal outlet pressure
Q1 = Flow

1/2 Symbole (according to NF EN 980)

 = Attention, instruction data sheet attached to be carefully read


2. AREA OF APPLICATION

The pressure-regulator **NOVA 40** is designed to be used for compressed medical gases except instrument gases. It is equipped at the outlet with a standardized medical terminal unit. It fits onto the valve of a compressed gas cylinder at 200 bar maximum. Its connection may be realized either by a nut and a stem on a cylinder valve or by a yoke for a pin-index valve. Intended for emergency medicine (ventilation) and oxygenotherapy.

3. FEATURES

3/1 Operating features

- P1 : 200 bar
- P2 : preset 4 bar in the standard version
- Flow : 50 l/min (at the outlet of the terminal unit)
- Opening of the relief valve : 7 bar
- Operating and warehousing temperature : - 20°C to + 60°C

3/2 Dimensional features

- Standardized inlets for cylinder, threads or yoke
- Rear inlet
- Weight : 940 g

The present data sheet only describes the standard version and its major options. It does not take into account the different versions that may be supplied.

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4. SAFETY

This equipment works with gas under pressure - ABSOLUTELY FOLLOW THESE INSTRUCTIONS.

- *Make sure that the cylinders are stored in a well ventilated area and distant from any source of heat and easily flammable materials.*
- *Do not smoke nearby.*
- *Make sure that the gas cylinders are correctly secured on a flat and horizontal surface, use a cylinder holder to limit the risks of falling.*
- *The relief valve supplied with and assembled on this equipment is exclusively intended for the protection of said equipment. A careful risk analysis regarding potential gas returns as well as pressures and flows of the downstream pipeline system must be carried out in order to determine the need of additional protections.*
- *If necessary, provide additional safety devices to prevent risks from overpressures or gas return.*
- *This equipment is provided with a relief valve.*
- *Slowly and progressively open the gas cylinder valves as well as all kinds of valves, (for oxygen, it will thus prevent the effects of overheating by compression).*
- *Protect the cleanliness of this equipment. Penetrations of particles during handling are responsible for most breakdowns.*
- *Do not sterilize the unit.*
- *Only use equipment in perfect condition.*
- *Only use equipment compatible with the intended gases, the required pressures and flow rates.*
- *Check and maintain a perfect gastightness of the equipment and pipes.*
- *Only use connections delivered with the equipment, do not use intermediate connections.*
- *Never work on an equipment or pipe under pressure or modify the setting of the safety relief valve.*
- *To maintain the equipment in good working conditions and in case of abnormal variations of performances, consult the paragraph "Maintenance".*
- *The marking "HP" on certain orifices indicates the presence of the upstream pressure (high pressure).*
- *Installations comprising high flow pressure regulation of CO₂ and/or N₂O may show some "icing". Consult your fitter about the necessity to install a heater.*
- *An unit that has been subjected to high temperature (> 100°C) must be discarded.*
- ***Never lubricate an equipment used with compressed gases and particularly oxidant gases (O₂, N₂O, NO...).***
- *So as to prevent any overpressure into the terminal unit, it is FORBIDDEN to alter the adjustment of the relief valve except during maintenance operations.*
- *In case of fall or chocs on the device, systematically recheck its functions and refer to the "Maintenance" chapter if necessary.*
- *Systematically protect the inlet and outlet of any non connected equipment.*

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
- Although presenting no particular difficulties, the commissioning of this device should only be handled by authorized and properly trained persons.

- As there is no alarm device on the equipment, the patient must be watched over.

- The equipment must be out of reach from children or irresponsible persons.

Not respecting these rules may cause severe damages.

5. INSTALLATION

5/1 Precautions to be taken before installation

- Verify that the packaging is not damaged.
- Upon delivery, check for any damage that may have occurred during transport. Verify that the content corresponds to your order and needs.
- The equipment delivered is free of grease. Care should be taken to prevent pollution during installation

5/2 Installation on the cylinder valve

During the first installation and during each cylinder change.

- Close the cylinder valve.
- Purge and then, if necessary, unconnect the downstream line.
- In the event of a cylinder change, unmount the **NOVA 40** by unscrewing from the cylinder's valve, the inlet connection nut or the yoke.
- Check the presence of the seal at the end of the inlet connection or on the yoke as well as the cleanliness of all the connections.
- Check that the inlet connection or the yoke match the outlet connection of the cylinder.
- Fully screw in the inlet connection nut or the yoke on to the cylinder valve.
- Connect the downstream line.
- If you have ordered a version without inlet connection, it is mandatory to use a standardized connector according to your local regulations. In this case, carefully follow the following steps :
 - Check the presence of the inlet filter into the body,
 - Screw the standardized connection corresponding to the gas into the body (torque : 65 N.m),
 - Check the airtightness of the connection when putting into service (foam reagent).

NOTE : Handtighten the connections equipped with O-rings (ex: cylinder connection), connections with other types of seals are to be fitted by means of a wrench (ex : pipe connection).

For safety reasons while using the equipment, oil or grease should be strictly prohibited while mounting. The symbol  and the indication "NO OIL" reminds the user of this precaution.

6. START-UP

6/1 Preliminary checks

- Always stand on the side and not in front of the **NOVA 40**.
- To avoid the effects of overheating by compression, open the cylinder valve slowly.

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- Carry on flushing so as to drain possible impurities.

- Check the absence of leaks with a foam reagent.

- Check the pressure of the gas cylinder on the high pressure gauge.

6/2 Terminal unit

- Connect the medical device onto the terminal unit of the **NOVA 40** by pushing, turning or screwing the probe of the device.
- Locking and airtightening are obtained according to your regional standard.
- The flow is immediately available when the probe is in place.

IMPORTANT : In absence of use for a long period, shut-off the line upstream. In any case, the **NOVA 40** should never be used as a shut-off valve.

6/3 Desinstallation from the cylinder

- To desinstall the **NOVA 40**, after closing the upstream circuit, purge the downstream part (through the terminal unit if need be).
- Do not unscrew or tighten the inlet connector, when the gauge is under pressure.
- Do not forget to disconnect the devices that may be attached to the **NOVA 40**.

7. MAINTENANCE

The safe use of this product requires the following regular precautions :

- Daily check of the good overall aspect and absence of leaks.
- Regular maintenance of this product according to its degree of use.

Servicing time every	5 years	3 years	1 year
Use	low or shelf	8h/day	24 h/day

- Discarding after 15 years in working or in storage or 5 repairs, whichever the shortest period.
- Maintenance kit : YM 1352

In no case should the maintenance periodicity of the regulator be over 5 years, and 1 year for the control of the flow through the terminal unit.

Servicing needs no particular tool but should be carried out only by authorized personnel.

Use only original spare parts. In case of important repairs (accidental damage, abnormal variations in performance), the equipment should be returned to the manufacturer.

The terminal unit does not require any particular maintenance outside a regular cleaning.

NT 0647 +