



**PRODUCT
INFORMATION
AND INSTRUCTIONS**

OM-700

OM-700S

LOTUS[®]
electronic oxygen conserver



TABLE OF CONTENTS

General Information	2
Important Safety Rules and Precautions	3-4
Introduction	5
Description of Parts and Controls	6-8
Assembly and Use	
<i>Non-Portable Use</i>	9
<i>Installing Alkaline Batteries</i>	10
<i>Installing the Battery Pack</i>	11
<i>Installing the System</i>	12
<i>Charging the Battery Pack (rechargeable NiMH Battery Pack)</i>	13
<i>Charging the Battery Pack (full charge)</i>	14-15
<i>Operating Instructions</i>	16
<i>Alarm Operation</i>	17
Oxygen Cylinder Duration	18
Care and Maintenance	
Disposal	19
Accessories	20
Troubleshooting Guide	21-22
Classifications	23
Specifications	24-25
Symbols Key	25
Limited Warranty	26
Important Information to Record	27
Notes	28



GENERAL INFORMATION

This manual provides information necessary to operate the LOTUS® electronic oxygen conserver with built-in regulator in accordance with a physician's prescription.

The LOTUS® conserver can be used with any CGA 870 post-valve cylinder (see Fig. A) at home or away from home to provide your specific oxygen requirements. It requires two (2) 1.5 volt AA alkaline batteries for operation, or, it can be used with a rechargeable battery pack that may be purchased separately.

Statements in this manual preceded by the following words are of special significance:

WARNING

Indicates there is a possibility of injury to you or others.

CAUTION

Indicates there is a possibility of damage to the device or to other property.

NOTE

Indicates points of particular interest or emphasis that allow for more efficient and convenient operation of the equipment.





FIGURE A
CGA 870
Post-Valve Cylinder

IMPORTANT SAFETY RULES & PRECAUTIONS

WARNING:

- Read and understand this manual before operating your LOTUS® electronic oxygen conserver.
- This device is not intended for use by patients who:
 - Breathe more than 40 breaths per minute,
 - Consistently fail to trigger equipment (i.e., mouth breathe).
-  **Smoking near oxygen equipment is strictly prohibited.** Keep cigarettes, matches, burning tobacco, and open flames, such as lighted candles, away from the area where the system is being stored or operated.
- Avoid creation of any spark, such as static electricity caused by any type of friction, near the oxygen equipment.

 **NOTE: Oxygen will not burn; however, it does vigorously accelerate the burning of any flammable material.**

- Never use oil, grease, or petroleum-based products on or near the system. Please wash and dry your hands properly prior to operating your oxygen equipment.
- Never use aerosol sprays near the oxygen equipment.
- Do not use in the presence of flammable anesthetic mixture.
- Keep your oxygen equipment at least 5 ft. (1.5 m) away from any electrical appliance.
- Be sure to turn off the oxygen supply by closing the cylinder valve when not in use.
- Do not use cannula tubing that is longer than 7 ft. (2.13 m).
- Do not use mask, pediatric, or other low-flow cannula tubing when operating the unit.



IMPORTANT SAFETY RULES & PRECAUTIONS

CAUTION:

- Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.
- Prevent water or other liquid substances from entering the unit.
- Prevent dust or any small particles from entering the unit.
- Do not expose the unit to extreme temperatures.
- Be sure to carry extra AA-size **alkaline** batteries in the event they are needed.
- Always maintain a backup supply of oxygen.
- Do not use humidifier bottles with the unit.
- Do not use if leaking or damaged.
- Refer repairs to authorized service personnel.


 **NOTE: Oxygen supplied by this equipment is supplemental only and is not intended for life support applications.**

Please contact your Home Care Provider
if you have any questions.

INTRODUCTION

The LOTUS® electronic oxygen conservers include a combination of a low-pressure regulator and an oxygen conserver. They are designed for use with ambulatory systems and are capable of delivering a precise amount of supplemental oxygen at the optimal point in the breathing cycle. Operationally, the LOTUS® conservers greatly increase efficiency in the delivery of oxygen, maximizing the beneficial effects and eliminating unnecessary oxygen waste.

When we breathe, approximately one-third of the time is spent inhaling and two-thirds exhaling. As a result, oxygen delivered by continuous flow is wasted during exhalation. By eliminating oxygen flow during exhalation, a two-thirds savings is possible. Additionally, the oxygen available during the very first part of inhalation contributes most to meeting oxygen needs. The LOTUS® conservers take advantage of these facts to provide maximum efficiency in the delivery of oxygen. These devices are designed to be an integral component of a lightweight, long-lasting ambulatory oxygen system.



DESCRIPTION OF PARTS & CONTROLS




FIGURE B
LOTUS® Conservor
Top View







FIGURE C
LOTUS® Conservor
Bottom View

DESCRIPTION OF PARTS & CONTROLS

- **Battery Test Indicator Light:** The LOTUS® conservers incorporate a visual indicator light (multi-color LED) that monitors the battery energy level inside the unit by displaying three colors.
 - 1) **Green** - adequate battery energy level for operation.
 - 2) **Amber** - lower battery energy level, but still in operational range.
 - 3) **Red** - low battery energy level. Two (2) spare AA-size alkaline batteries, an extra, fully-charged battery pack or access to a recharger should be readily available.
 - 4) **Flashing Red** - Replace batteries immediately.
- **Cylinder Adjustment Knob:** This is used to attach the unit to any CGA 870 post-valve cylinder.
- **Selector Knob:** This enables the user to select the prescribed oxygen flow setting. It also enables the user to monitor the battery energy level. When not in use, the switch should be turned to the “OFF” position.
- **Battery Holder:** This compartment holds two (2) AA-size alkaline batteries or one (1) rechargeable NiMH battery pack (RB-100).
- **Oxygen Pressure Gauge:** This enables the user to monitor the contents of the compressed oxygen cylinder.
- **Oxygen Supply Outlet:** Use this fitting to attach a standard cannula.
- **Cylinder Alignment Pins:** When assembling the unit, these parallel pins must go in the holes on the post valve.
- **CGA Seal (*Gasket FR-870G, brass and viton*) or equivalent:** This creates the interface between the valve and the LOTUS® conserver. Besides offering a rugged interface, it also surrounds the oxygen path in a ring of stainless steel or brass.
- **Oxygen Inlet Connection:** Interface between cylinder valve and conserver that allows oxygen to flow into the regulator.
- **Battery Charger Connector Port:** This port allows the user to charge the rechargeable NiMH battery pack from a standard AC outlet. (Requires the use of an optional manufacturer-specified rechargeable battery pack, charger, and adapter. See the “Accessories” section on page 20 for details.)
- **Charger Indicator:** This visual indicator light (green LED) indicates when the rechargeable battery pack is being charged. This light will also come on if the selector knob is set to battery position () when using the NiMH battery pack.

WARNING: Use only a manufacturer-specified gasket. Other gaskets may not be oxygen compatible and may cause an oxygen leak, creating an increased fire risk.



DESCRIPTION OF PARTS & CONTROLS

- **Bypass Switch:** This is a pneumatic switch that enables the user to switch from pulse mode (oxygen delivery on demand) to continuous flow mode in the unlikely event of unit malfunction. Continuous flow mode is not powered by the battery and can be utilized without one.

☞ **NOTE:** Remember that in continuous flow mode, the oxygen will be consumed at a much faster rate. Return to another source before depleting the oxygen cylinder.



FIGURE D
Close-up of Bypass Switch

CAUTION: The Bypass Switch is designed for emergency use only. In the event that it is necessary to operate the unit in this continuous flow mode, **DO NOT** obstruct the flow of oxygen from the Oxygen Supply Outlet by placing your finger over the outlet or blocking the flow through the oxygen tubing in any way. Doing so may render the unit inoperable and/or damage the sensor in the unit.

☞ **NOTE:** The continuous flow function on the LOTUS® conserver is factory preset at 2 lpm. It can be adjusted between 0.5 and 6 lpm by your Home Care Provider.


ASSEMBLY AND USE

- Make certain that your hands are free of oil, grease, and other contaminants.
- Inspect the unit to insure that it has a CGA 870 brass and viton (or equivalent) gasket in good working condition attached to the inlet nozzle.
- Secure the cylinder in an upright position.
- Inspect the post valve of the cylinder and the LOTUS® conserver to ensure they are free of contaminants. If any indication of damage or contamination is detected, DO NOT use the equipment and contact your Home Care Provider.

WARNING: Use ONLY a manufacturer-specified gasket. Other gaskets may not be oxygen compatible and may cause an oxygen leak, creating an increased fire risk.

NON-PORTABLE USE:

The LOTUS® conservers are designed to extend the life of portable oxygen supplies when away from the primary source. While the LOTUS® conservers may be used with stationary oxygen cylinders they should be used only while awake and reasonably attentive.



ASSEMBLY AND USE

INSTALLING ALKALINE BATTERIES: (Standard AA batteries)

- STEP 1:** While holding the unit with one hand, press in and lift up on the battery holder door latch.
- STEP 2:** Press the batteries into the slot, making sure they are placed over the battery release strap and inserted in the proper direction [Fig. E]. They should fit snugly into the unit. Be sure that the batteries make firm contact at the positive (+) side.






FIGURE E
Installing AA alkaline
batteries

- NOTE:** The LOTUS® conserver is packaged with batteries inside. A special seal is used at the positive terminal to prevent battery oxidation. If it has not already been removed by your Home Care Provider, be sure to remove the seals before using the unit for the first time.
- STEP 3:** Close the battery holder cover by snapping it back into place.
- STEP 4:** Verify the battery level by turning the selector knob to the battery position ().
- NOTE:** The battery indicator light is inoperative except when the selector knob is in the battery position (); however, the red light will flash in any setting if the batteries need to be replaced. Always check the battery energy level before using. With normal use of eight (8) hours per day at setting two (2), a battery should last approximately six (6) months.

ASSEMBLY AND USE

INSTALLING THE BATTERY PACK:

(Rechargeable NiMH AA Battery Pack - Product Code RB-100)

STEP 1: While holding the unit with one hand, press in and lift up on the battery holder door latch.

STEP 2: Press the battery pack into the slot so that you can read the words “This Side Up” on the battery label and so the arrow is pointing toward the battery door [Fig. F]. The battery pack should fit snugly into the unit.





FIGURE F
Installing the
NiMH AA
battery pack.

STEP 3: Close the battery holder cover by snapping it back into place.

STEP 4: Verify that the battery is charged by turning the selector knob to the battery position ().

If this is your first time using the battery pack, you must first charge it for approximately 4-6 hours before using.

NOTE: With normal use of eight (8) hours per day at setting two (2), a battery charge should last approximately one (1) month. Because of the nature of NiMH batteries, the battery pack will discharge in 2-3 months if not in use. If this occurs, follow the recharging procedure before using.

ASSEMBLY AND USE

INSTALLING THE SYSTEM:

- STEP 1:** Loosen the cylinder attachment knob.
- STEP 2:** Lower the LOTUS® conserver over any CGA 870 post-valve cylinder with the alignment pins toward the holes on the cylinder neck [Fig. G].
- STEP 3:** Align the two (2) pins and gasket with the corresponding holes on the cylinder post valve.
- STEP 4:** While holding the unit in place, tighten the cylinder adjustment knob by turning clockwise [Fig. G].

 **NOTE:** Tighten only by hand. The use of a tool to tighten the knob may damage the unit.

CAUTION: If you are unable to eliminate leaks by manually tightening the cylinder adjustment knob, replace the CGA 870 gasket. If leaks persist, the unit must be returned for service.

- STEP 5:** Attach a standard cannula to the oxygen supply outlet.



FIGURE G
Attaching the LOTUS® conserver to the cylinder

ASSEMBLY AND USE

CHARGING THE BATTERY PACK (rechargeable NiMH battery pack):

The LOTUS® recharger kit (Product Code RK-100) is comprised of three (3) components [Fig. H]:

- 1) Rechargeable NiMH AA Battery Pack (RB-100)
- 2) Recharger (RC-100)
- 3) Recharger Plug Adapter, U.S. 110VAC (RA-110)




FIGURE H

ASSEMBLY AND USE

CHARGING THE BATTERY PACK (full charge):

STEP 1: Connect the Plug Adapter to the Recharger input socket [Fig. H1]. Plug the Adapter into the AC power outlet. The red indicator light will flash three (3) times.

STEP 2: Connect the Recharger output connector into the LOTUS® Battery Charger Connector Port [Fig. H2]. The battery pack must be installed in the LOTUS® conserver before charging. If the recharger components are connected properly, both the red light on the recharger and the green charger indicator light on the LOTUS® conserver will illuminate.

STEP 3: The battery is now charging. The initial charge may take up to six (6) hours. For all future charges, the charge time will be approximately four (4) hours. When the battery is fully charged, the red light on the recharger will turn off and the green indicator light will turn on.

STEP 4: Remove the recharger from the LOTUS® Battery Charger Connector Port. The LOTUS® conserver is now ready to be used as part of your ambulatory oxygen system.




FIGURE H1

Connect Plug Adapter to the Recharger Input Socket





FIGURE H2

Connect the Recharger Output into the LOTUS® Battery Charger Connector Port

ASSEMBLY AND USE

CHARGING THE BATTERY PACK (full charge): continued

- If the recharger is removed from the power source, you must wait fifteen (15) seconds for the recharger to reset before reconnecting it to the LOTUS® conserver and beginning the charging process. Also, leaving the recharger connected to the LOTUS® conserver after the green light on the recharger illuminates (indicating charge is complete) will not damage the system.
- After the initial charge, you may use the LOTUS® conserver while charging. After the battery module is fully charged, the recharger will return to trickle charge mode to top off the battery as it is being depleted.
- For full instructions, specifications, and warnings, please refer to the instruction sheet included with the recharger kit.





ASSEMBLY AND USE

OPERATING INSTRUCTIONS:

STEP 1: Make sure the Bypass Switch is set to Pulse Mode, and the LOTUS® conserver is set to the “OFF” position.

STEP 2: To reduce the risk of rapid oxygen recompression and fire, **OPEN THE CYLINDER VALVE SLOWLY** and completely so that the pressure gauge moves slowly as it indicates the cylinder pressure.

STEP 3: Listen for leaks. If a leak is present, close the cylinder valve, check the CGA seal, and reinstall. If the leak persists, **DO NOT USE THE EQUIPMENT**. Contact your supplier for repair.

STEP 4: Turn the setting selector knob to the battery position (⊕⊖). Check the battery indicator:

- 1) **Green** - Good
- 2) **Amber** - Lower battery energy level, but still in the operational range.
- 3) **Red, Flashing Red, or No Indicator** - Replace the battery.

NOTE: The battery indicator light will flash amber for one (1) second and the audible alarm will beep once before displaying the battery energy level.

STEP 5: Select the setting on the LOTUS® conserver (1-6) that corresponds to the appropriate delivery setting.

STEP 6: Place the nasal cannula into position with the prongs in the nostrils and begin breathing.

The LOTUS® conserver will now start to deliver oxygen. The amount of oxygen delivered is determined by the setting. A “click” sound may be heard each time the unit delivers a pulse of oxygen. Adequate oxygen delivery will be achieved because of the precise time in the breathing cycle that the pulse of oxygen is delivered.

NOTE: To help prevent possible damage to the unit, keep the LOTUS® conserver in a carrying bag. Several bags are available for use with different cylinder sizes and configurations.

STEP 7: When finished using the system, turn off the oxygen supply cylinder valve and continue breathing through the nasal cannula until no further oxygen is detected.

STEP 8: Remove the nasal cannula and turn the selector switch to the “OFF” position.

STEP 9: When not in use, store in a clean, dry location.

ASSEMBLY AND USE

ALARM OPERATION:

LOW BATTERY ALARM (ALL MODELS)

When the LOTUS® conserver senses a low battery level while the device is being operated in a delivery setting (between 1 and 6), the following will occur:


- The battery indicator LED will flash red.
- An audible alarm will beep every fifteen (15) seconds until the battery is replaced, recharged, or depleted completely.

👉 **NOTE:** If the battery level drops while the device is in the battery position (←→), the red indicator LED may not flash.

BREATH SENSOR ALARM (Model OM-700S ONLY)

This alarm monitors your breathing. If no breath is sensed for 40 seconds, the following will occur:

- The battery indicator LED will flash red.
- An audible alarm will beep every second until a breath is sensed.



OXYGEN CYLINDER DURATION

Because the total delivery of oxygen via the LOTUS® conserver is related to breathing rates, it is user adaptive in that the total oxygen delivered per minute will automatically adjust with user need, as expressed by increased or decreased breathing rates. For example, at all settings, twice as much oxygen per minute will be delivered if one breathes twenty (20) times per minute as compared with ten (10) times per minute. The table below provides useful information to be used as a guide.

	SETTING	1	1.5	2	2.5	3	4	5	6	Continuous Flow 2lpm
CYLINDER TYPE	CYLINDER VOLUME	Estimated Cylinder Duration in Hours (Based on 20 breaths per minute)								
M2	36 liters	3.0	2.0	1.5	1.2	1.0	0.8	0.6	0.5	0.3
M4(A)	113 liters	9.4	6.3	4.7	3.8	3.1	2.4	1.9	1.6	0.9
M6(B)	164 liters	13.7	9.1	6.8	5.5	4.6	3.4	2.7	2.3	1.4
ML6	171 liters	14.3	9.5	7.1	5.7	4.8	3.6	2.9	2.4	1.4
M9(C)	246 liters	20.5	13.7	10.3	8.2	6.8	5.1	4.1	3.4	2.1
D	425 liters	35.4	23.6	17.7	14.2	11.8	8.9	7.1	5.9	3.5
E	680 liters	56.7	37.8	28.3	22.7	18.9	14.2	11.3	9.4	5.7

CARE AND MAINTENANCE

The LOTUS® conserver is designed for a long and accurate life; however, as with any electronic device, prudent care is required. The unit should be kept clean and free from moisture and dust, as well as extreme temperature. Do not expose the unit to water, such as when bathing or swimming. It is advisable to keep the device in a carrying bag to afford a degree of protection. Clean the outside cover periodically by wiping with a lint-free cloth. Pay special attention to the oxygen inlet and outlet to make sure they remain free of dust, etc. If the oxygen inlet connection becomes contaminated with dirt, oil or grease, **DO NOT USE OR ATTEMPT TO CLEAN**. Contact your supplier for service or repair.

Cannula tubing is a disposable accessory that should be replaced periodically following normal usage. Disposable tubing should be disposed of in accordance with local ordinances and/or regulations for disposal. Replacements are available through your Home Care Provider (OC-401, case of (50) 4 ft [1.22 m] cannulas, or equivalent).

DISPOSAL:



The Unit: Do **not** dispose of the unit in the household waste. Consult an authorized electronic waste recycling company for the proper disposal of the unit. You can find out their address from your environmental officer or from your local council.



Disposal of Batteries/Battery Pack: Do **not** dispose of used batteries/battery packs in the household waste. Contact your distributor or a public waste disposal authority.



ACCESSORIES

RB-100 Rechargeable NiMH AA Battery Pack

RC-100 Recharger

RA-110 Recharger Plug Adapter, U.S., 110VAC

RK-100 Recharger Kit, U.S. Configuration (includes battery pack, recharger, and plug adapter)

 **NOTE: Only manufacturer-specified battery packs and rechargers may be used with the LOTUS® conserver. These accessories are available from your Home Care Provider.**

TROUBLESHOOTING GUIDE

PROBLEM	PROBABLE CAUSE	SOLUTION
Unit does not pulse.	Dead batteries.	Replace or recharge batteries.
	Batteries installed incorrectly (reversed).	Make sure battery polarity is correct.
	Dirty battery holder contacts.	Remove the batteries. Use rubbing alcohol and a cotton swab to clean contacts.
	Cylinder valve is closed.	Turn the cylinder valve to the "ON" position.
	Cylinder is empty.	Check the oxygen gauge. Replace the cylinder, if empty.
	Oxygen cannula is blocked or kinked.	Remove kinks. Clean or replace, if necessary.
Short battery life.	Non-alkaline batteries are used.	Make sure the batteries inside the unit are alkaline.
	Battery pack was not fully charged before it was inserted or before the device was used.	Leave recharger adapter plugged into the connector port until the red light on the recharger turns off and the green light turns on.
Red indicator on recharger does not blink three (3) times upon initial connection.	Plug adapter is not properly connected to the electrical wall outlet.	Make proper connection.
	Plug adapter is not properly connected to the recharger input socket.	Make proper connection.
	Electrical wall outlet is not receiving power.	Check power source.
	Faulty plug adapter/recharger.	Discontinue use and contact your Home Care Provider.
LOTUS® does not beep upon start-up.	Dead batteries.	Replace or recharge the batteries.

TROUBLESHOOTING GUIDE

PROBLEM	PROBABLE CAUSE	SOLUTION
The red indicator on the recharger does not illuminate to indicate that the battery pack is charging.	Recharger is not properly connected to the battery charger connector port.	Make proper connection.
	Battery pack is not installed correctly inside the battery holder compartment.	Make sure that the battery pack is installed correctly. Refer to “Installing the Battery Pack” section of this book on page 11.
	Dirty battery holder contacts.	Remove the battery pack. Use rubbing alcohol and a cotton swab to clean all contacts.
	Battery strap is caught between the battery pack and the battery holder side contacts.	Make sure that the battery strap is not caught between the battery pack and the battery holder side contacts.
	A non-approved battery pack or alkaline batteries are installed.	Install the proper battery module (RB-100).
	Defective battery pack.	Discontinue use and contact your Home Care Provider.
Breath sensor alarm does not operate or operates intermittently.	Non-alkaline batteries are installed.	Install alkaline batteries.
	Batteries are low/dead.	Install new batteries.
	Batteries are non-functioning.*	Install new batteries.

*Because of variations in Duracell®’s manufacturing process, the alarm on the LOTUS® Model OM-700S may not work properly when using Duracell® batteries.

 **NOTE: For further troubleshooting tips on the rechargeable battery module or recharger, see the instruction sheet packaged with the recharger kit.**


Non-functioning units are subject to warranty provisions and the manufacturer repair/return policy. If necessary, call your Home Care Provider.

 **NOTE: Do not attempt to open the electronic compartment of the unit. If the case is opened or tampered with, the warranty is void.**

CLASSIFICATIONS

The LOTUS® conserver is classified as:

- Internally Powered Equipment - EN 60601-1
- Ordinary Equipment - Enclosed equipment without protection against ingress of water, EN 60601-1
- Not suitable for use in the presence of flammable anesthetic mixture with air, oxygen or nitrous oxide
- Complies to EN 60601-1-2 electromagnetic testing
- Type B Equipment
- Class II, per FDA 21 CFR Part 868.5905





SPECIFICATIONS


Oxygen Delivery:	<u>Switch Position</u>	<u>Liter Flow Equivalency</u>
	1	1
	1.5	1.5
	2	2
	2.5	2.5
	3	3
	4	4
	5	5
	6	6
Continuous Flow Emergency Bypass Setting:	Factory preset at 2.0 ± 0.5 lpm, provider adjustable between .5 and 6 lpm	
Required Operating Pressure:	200 psi to 3000 psi (13.8 bar to 206.8 bar)	
Regulator:	Built-in, 25 ± 5 psi ($1.7 \pm .3$ bar)	
Voltage:	2.4 VDC to 3.0 VDC	
Battery:	(2) x 1.5 VDC Alkaline AA-sized batteries or (1) rechargeable NiMH battery pack	
Battery Indicator:	Green:	Good
	Amber:	Lower battery energy level, but still in the operational range
	Red:	Replace
	Flashing Red:	Replace immediately
Power Consumption:	<10 mA at idle time and maximum of 600 mA during delivery time	
Dimensions:	Approximately 7" L x 2.5" H x 4" W (17.8 cm L x 6.4 cm H x 10.2 cm W)	
Weight:	Approximately 14.5 oz. (408 grams)	

SPECIFICATIONS

Operating Temperature:	32°F to 122°F (0°C to 50°C)
Operating Relative Humidity:	15% to 95%
Operating Altitude:	0 to 10,000 feet (0 to 3,048 meters)
Storage/Transportation:	Cold temperature: Maximum -40°F (-40°C), 1% RH Hot temperature: Maximum 145°F (63°C), 44% RH
Shock:	Not to exceed IEC 601-1 requirements
Vibration:	Not to exceed IEC 68-2-6, IEC 68-2-34


SYMBOLS KEY:

 : No smoking or open flames

 : Consult accompanying documents

| : System ON

○ : System OFF

 : Type B equipment





LIMITED WARRANTY

The LOTUS® oxygen conserver has been carefully manufactured and inspected and is warranted to be free from defects in workmanship and materials. Under this warranty, CHAD Therapeutics' obligation shall be limited to the replacement or repair of any such units or parts that prove, by CHAD's inspection, to be defective within two years from the date of purchase. Any abuse, operation other than the intended use of the product, negligence, accident or repair by other than authorized service professionals shall immediately void this warranty. This warranty does not extend to the cannula or battery.

CHAD Therapeutics will not accept damages or charges for labor, parts or expenses incurred in making field repairs, except upon written authorization prior to such action.

The foregoing warranty is exclusive and in lieu of all other express warranties. Implied warranties, if any, including but not limited to the implied warranties of merchantability and fitness for a particular purpose, shall not extend beyond the duration of the express warranty provided herein. In no event shall CHAD Therapeutics be liable for loss of use or profit or other collateral, special or consequential damages.

IMPORTANT INFORMATION TO RECORD

Your Name: _____

Date You Received Your Unit: _____

Prescribed Oxygen Flow Setting:

- At Rest: _____
- During Exercise: _____


Home Care Provider's Name: _____


Home Care Provider's Phone Number: (_____) _____

Physician's Name: _____


Physician's Phone Number: (_____) _____

Notes: _____





NOTES



Manufactured by Inovo, Inc.

2975 Horseshoe Drive South, Suite 600
Naples, FL 34101

Toll-free: 800-423-8870

Local: 239-687-1280

Fax: 239-687-1285

International Phone: +01-818-701-1008

International Fax: +01-818-701-1913

www.chadtherapeutics.com

Authorized Representative
According to MDD 93/42/EEC
MDSS

Schiffgraben 41
30175 Hannover, Germany

United States Patent Number # 5,134,886
PM-00010/09/K

Printed in the U.S.A.

LOTUS and CHAD are registered trademarks of Inovo, Inc.



ISO 13485
Registered