

## 'B' Size Cylinder System



## A Truly Portable Medical Oxygen System

The 'B' size cylinder is the smallest and lightest medical oxygen cylinder in Australia. This truly portable oxygen cylinder ends the limitation of mobility and enables you to be free to travel without boundaries.



**Truly Portable**



**Lightweight**



**Airline Approved**



**Extended Usage**



**Economical**



## What is the 'B' size cylinder system?

The 'B' size cylinder system consists of a 160L medical oxygen cylinder and an oxygen conserving device (OCD) contained within a robust carry bag. This lightweight portable medical oxygen system is designed to enhance your quality of life, delivering up to seven<sup>1</sup> hours of oxygen at 2 LPM when used with an OCD.

## What is an OCD?

An OCD is a special battery-operated demand regulator, releasing oxygen when you inhale and not when you exhale. This 'pulse dose' technology extends the time that an oxygen cylinder will last.

## Why should I use the 'B' size cylinder Portable Medical Oxygen System?

### Dual flow settings

- ✓ Two flow settings: Continuous or Pulse
- ✓ Pulse settings from 1 to 6 LPM

### Extended Usage

Estimated oxygen duration<sup>1</sup> in hours for a 'B' size cylinder (160L) using an AirSep® Impulse™ Elite OCD<sup>2</sup>

|               | 1 LPM   | 2 LPM   | 3 LPM   |
|---------------|---------|---------|---------|
| <b>Mode A</b> | 15 hrs  | 7.6 hrs | 5 hrs   |
| <b>Mode B</b> | 8.1 hrs | 4.1 hrs | 2.7 hrs |
|               | 4 LPM   | 5 LPM   | 6 LPM   |
| <b>Mode A</b> | 3.8 hrs | 3 hrs   | 2.5 hrs |
| <b>Mode B</b> | 2 hrs   | 1.6 hrs | 1.6 hrs |

### Lightweight

- ✓ Full cylinder with OCD & carry bag weighs only 2kgs

### Economical

- ✓ OCD operates on 1.5V alkaline battery
- ✓ Battery lasts up to 12 months on normal usage

### Ease of Travel

- ✓ Convenient to use and transport in private vehicle
- ✓ Domestic airline approved<sup>3</sup>



| Specifications           |                                               |
|--------------------------|-----------------------------------------------|
| <b>'B' size cylinder</b> |                                               |
| <b>Weight</b>            | 1.3kgs full without OCD<br>2kgs full with OCD |
| <b>Dimensions</b>        | 39cm H x 8.2cm Diameter                       |

## Oxygen without Boundaries

Oxygen without Boundaries is a program designed to help you choose a portable medical oxygen therapy solution to suit your lifestyle.

With great benefits, our rental program can assist you in your personal mobility requirements. The 'B' size cylinder system is available for rental on a monthly basis (OCD also available for purchase).

<sup>1</sup> Based on 20 breaths per minute; <sup>2</sup> Duration of use will differ slightly when a different model of OCD is used

<sup>3</sup> Subject to conditions imposed by the carrier in its Travel Clearance Guideline

Contact us now to find out more about our services on 1300 36 02 02.

**NSW State Office**  
Unit 5/476 Gardeners Road  
ALEXANDRIA NSW 2015  
Tel: 02 9364 7474  
Fax: 02 8338 9797

**QLD State Office**  
25 Campbell Street  
BOWEN HILLS QLD 4006  
Tel: 07 3252 8039  
Fax: 07 3252 4850

**VIC & TAS State Office**  
40 Bunnet Street  
NORTH SUNSHINE VIC 3020  
Tel: 03 9310 1200  
Fax: 02 9364 7482

**SA State Office**  
95 Fullarton Road  
KENT TOWN SA 5067  
Tel: 08 8331 1648  
Fax: 08 8331 1849

