

AIR LIQUIDE *Healthcare*:

> Over 20 years of respiratory care experience worldwide

> 40 locations around Australia

Portable Medical Oxygen Rental Program

Oxygen Without Boundaries

Equipment Information

The 'B' size cylinder with OCD



Exclusive to AIR LIQUIDE *Healthcare*, this portable cylinder system is suitable for those of you who require mobility for occasional and short outings.

Specifications

'B' size cylinder	
Weight	<2kg with full cylinder and an OCD
Dimensions	39cm H x 8.2cm Diameter

PHILIPS RESPIRONICS EverGo™ Portable Concentrator



Suitable for those of you who require a portable unit that allows you to travel or work further and longer without feeling short of breath. The system comes with a trolley for increased transportability.

Specifications

Flow Setting	1-6lpm pulse flow
Weight	4.8kg with 2 batteries excludes trolley
Dimensions	16.8cm H x 27cm W x 10.7cm D
Battery Duration	up to 12hrs with 2 batteries

AirSep FreeStyle™ Portable Concentrator



Ideal for those of you who need a truly portable medical oxygen system to enjoy an active lifestyle. The unit's small size and light weight (2kgs) allows you to wear it while on the move!

Specifications

Flow Setting	1-3lpm pulse flow
Weight	Unit 2kgs / Air Belt™ 0.8kg optional
Dimensions	21.8cm H x 15.5cm W x 9.1cm D
Battery Duration	Internal - up to 3.5hrs with Air Belt™ - up to 10hrs

State Offices

New South Wales

State Office

ALEXANDRIA

Unit 5/476 Gardeners Rd
Alexandria 2015
Tel: 02 9364 7474

South Australia

State Office

KENT TOWN

95 Fullarton Rd
Kent Town 5067
Tel: 08 8331 1648

Queensland

State Office

BOWEN HILLS


25 Campbell St
Bowen Hills 4006
Tel: 07 3252 8039

Victoria & Tasmania

State Office

NORTH SUNSHINE

40 Bunnett St
North Sunshine 3020
Tel: 03 9310 1200



1300 36 02 02

www.airliquidehealthcare.com.au



This Paper is Carbon Neutral.



ALH 1ST ALH 2050-B 09/10



> Introducing Oxygen without Boundaries – A Truly Portable Oxygen Therapy Program

3 EASY STEPS TO START

Oxygen patients can now experience total freedom and reclaim their mobility thanks to the unique Portable Oxygen program from AIR LIQUIDE *Healthcare*. We empower you to get the solution you desire with a selection of equipment packages that can be rented or purchased at competitive rates. Experience the benefits of Oxygen without Boundaries today!

Truly Portable and Transportable Systems

Our portable medical oxygen solutions are selected to suit your individual mobility needs.

- ✓ For occasional outings
- ✓ To travel or work longer and further without feeling breathless
- ✓ To maximize your active lifestyle

Flexible Rental Options

Our rental program is designed to meet your personal mobility needs.

- ✓ Flexible terms: 1, 3 or 6 months
- ✓ Competitive rental rates
- ✓ Equipment trial available (up to 3 days)
- ✓ Purchase option available
- ✓ No deposit required

Peace of mind

Our large network of service outlets can provide comprehensive support in Australia.

- ✓ Therapy support
- ✓ Equipment repair services
- ✓ Loan units during equipment servicing
- ✓ Oxygen equipment assistance for travelling customers
- ✓ Equipment approved by a number of airlines

1. Review Prescription

Portable oxygen equipment uses “pulse dose technology”, and is suitable only for customers who can maintain adequate oxygen saturations on these devices. Please talk to your respiratory physician to find out if you are clinically suitable for this technology.

2. Start a Trial

Most portable medical oxygen systems operate on a pulse dose technology, which delivers medical oxygen only at the inspiratory stage of your breathing cycle. Since every person's breathing pattern is different, especially during exercising, it is therefore very important to trial the portable medical oxygen systems before making your rental or purchase decisions.

During our 3-day trial program you can:

- ✓ Verify with your physician if the pulse dose technology is suitable for you
- ✓ Determine the right level of oxygen flow settings for you at rest and during walking and exercising
- ✓ Decide if the chosen equipment is suitable for your lifestyle needs

Please note that we do not recommend the use of portable medical oxygen systems during sleep.

3. Choose Your Package

‘B’ Size Portable Medical Oxygen Cylinder System:

The ‘B’ size medical oxygen cylinder can be rented on a monthly basis with a standard regulator, or with an Oxygen Conserving Device (OCD), which can extend the duration of use by 3 or 6 times. A lightweight and comfortable carry bag with optional back strap can also be rented for your convenience.

Portable Medical Oxygen Concentrators:

AirSep® FreeStyle™ and Philips Respironics EverGo™

We have designed three rental packages to suit your individual mobility needs. You will also receive great benefits by choosing one of our rental packages. NO DEPOSIT required!

Rental Package Benefit	Gold*	Silver*	Bronze*
Minimum Rental Term	6 months	3 months	1 month
Education	✓	✓	✓
On-going Support	✓	✓	✓
Equipment Repair Services	✓	✓	✓
Replacement Unit during Repair	✓	✓	✓
Assistance in travelling with oxygen	✓	✓	✓
FREE starter kit	✓	✓	—
FREE Cannula (1 per month)	✓	✓	—
Homecare equipment discount	10% off RRP	5% off RRP	—

Our Portable Medical Oxygen Concentrators are also available for purchase.