

# Nitrous Oxide Sedation Equipment for Dental Practices

Specifically designed  
nitrous oxide sedation  
equipment and  
administration systems



AIR LIQUIDE *Healthcare* offers a range of equipment and administration systems specifically designed for dental surgery set ups.

- Range of options catering for a variety of needs
- Solutions compliant with industry recommendations

## NITROUS OXIDE / OXYGEN MIXERS

### ● Fraser MDM<sup>®</sup>

- Mechanical mixer
- Auto-compensation system
- Oxygen fail-safe and flush

Part # 583S



### ● Digital MDM<sup>®</sup>

- Digital mixer
- Audible and visual alarms
- Auto-compensation system
- Oxygen fail-safe and flush

Part # 3360S



## MOUNTING OPTIONS

### ● 2-cylinder mobile trolley

- Fully mobile stainless steel unit
- Accommodates 2 D-size cylinders
- Provides complete mobility
- Ideal for surgeries without piped gas



H: 88 cm  
W: 45 cm  
Depth: 30 cm

Part # 592S

### ● Telescopic pedestal

- Mobile stainless steel unit
- Adjustable height
- Mobile solution for surgeries equipped with piped gas



H (min): 70 cm  
H (max): 98 cm  
Radius: 65 cm

Part # 3335S

### ● Wall mounting

- 2 options available: 8" arm or telescopic 15" arm
- Discreet and aesthetic
- Can be stored in a cupboard
- Suitable for surgeries equipped with piped gas



8" arm

Part # 1268S



15" telescopic arm

Part # 605S

## PATIENT CIRCUITS & NASAL HOODS

### ● Matrix Scavenger

- Complete scavenging patient circuit
- Fully autoclavable
- Comes with a large nasal hood



Part # 3570S

### ● Reusable nasal hoods

- Suitable for Matrix Scavenger
- Fully autoclavable
- 3 sizes available



Small: Part # 848S  
Medium: Part # 847S  
Large: Part # 846S

### ● Disposable scented nasal hoods Dynomite<sup>™</sup>

- Suitable for Matrix Scavenger
- 5 fresh fragrances
- 3 sizes available



Small: Part # 1546S  
Medium: Part # 1547S  
Large: Part # 1548S

## OTHER EQUIPMENT & ACCESSORIES

### ● Breathing Bag

- 3 litre capacity
- Latex-free



Part # 852S

### ● Handwheel and hose assemblies



Oxygen  
Part # 3459SK



Nitrous oxide  
Part # 3460SK

### ● Regulators



Oxygen  
bull nose  
Part # 1801S



Nitrous oxide  
pin index  
Part # 944S

### ● 3" mounting pin

For use with 2 cylinder mobile trolleys and wall mounting options



Part # 1522S



For further information, please contact your AIR LIQUIDE Healthcare representative:  
Toll free 1300 36 02 02 / <http://www.airliquidehealthcare.com.au>