

**THE RESPIRONICS® FREEDOM SERIES™ —  
A NEW FAMILY OF PRODUCTS.  
AND A WHOLE NEW WAY TO LOOK  
AT OXYGEN THERAPY.**

The EverGo Portable Oxygen Concentrator is part of the Freedom Series of oxygen therapy products – an entire family of products designed to deliver more freedom – freedom from worry, freedom from discomfort, freedom from complicated maintenance. In short, the Freedom Series frees you from some of the more demanding aspects of oxygen therapy, so you can go on with your life.



**FEWER TIME CONSTRAINTS FOR EVEN MORE FREEDOM**

When you don't need to wait at home for oxygen deliveries, you are free to explore your world. Weighing under 10 pounds, and cleared for air travel by the FAA,\* EverGo is ideally suited for patients who are active and love to travel. That leaves you free to go wherever you want – whenever you want.

\* Note: Airline policies for traveling with an oxygen concentrator vary. Check with your airline before flying with EverGo.

**For more information about the EverGo™ Portable Oxygen Concentrator go to [evergo.respironics.com](http://evergo.respironics.com) or [www.evergo.respironics.eu](http://www.evergo.respironics.eu)**

**EVERGO™ PORTABLE  
OXYGEN CONCENTRATOR**

Basic system includes: EverGo Portable Oxygen Concentrator, Carrying Case (choice of color), Two Rechargeable Lithium Ion Batteries, DC and AC Power Cords, Accessory Case, and User Manual.

**PRODUCT**

EverGo System w/Black Carrying Case, International	900-020
EverGo System w/Navy Blue Carrying Case, International	900-022
EverGo System with Black Carrying Case, Canada	900-010
EverGo System with Blue Carrying Case, Canada	900-012

**ACCESSORIES**

EverGo Carrying Case – available in two colors	
<b>BLACK</b>	
Canada	900-107
International	900-123
Japan	900-125
<b>NAVY BLUE</b>	
Canada	900-108
International	900-124
Japan	900-126
Rechargeable Lithium Ion Battery	900-102
Accessory Case	900-250
Battery Charger	900-103
Gross Particle Filter	900-150
Mobile Cart	900-100
DC Power Cord	900-101
AC Power Supply	900-105
* Note: power cord alone is orderable based on your country	



Customer Service: +1-724-387-4000  
 Respironics Europe, Africa, Middle East: +33-1-47-52-30-00  
 Respironics Germany: +49-81-52-930-60  
 Respironics France: +33-2-51-89-60-00  
 Respironics Italy: +39-03-62-63-43-1  
 Respironics UK: +44-800-1300-845  
 Respironics Switzerland: +41-62-745-1750  
 Respironics Asia Pacific: +852-3194-2280  
 Respironics Australia: +61-2-9666-4444  
[www.evergo.respironics.com](http://www.evergo.respironics.com)

CAUTION: U.S. federal law restricts these devices to sale by or on the order of a physician.  
 Respironics is a registered trademark and EverGo is a trademark of Respironics, Inc. and its affiliates.  
 ©2007 Respironics, Inc. and its affiliates. All rights reserved.

Broudy CY 10/18/07 MCI 4101357 PN 1002329



**WE GAVE IT AN  
8-HOUR BATTERY LIFE —  
BECAUSE FREEDOM  
SHOULDN'T BE RUSHED.**

**EverGo™  
Portable Oxygen Concentrator**




# INTRODUCING THE EVERGO™ 8-HOUR, PORTABLE OXYGEN CONCENTRATOR

The EverGo Portable Oxygen Concentrator is the new oxygen therapy solution that truly gives the most active oxygen therapy patients the freedom to go wherever they want – whenever they want. Its 8-hour battery life\* and 1,050 ml/min oxygen capacity means that you can test your limits – without being limited by the device.

## LARGER OXYGEN CAPACITY LETS YOU STAY ACTIVE

Other portable concentrator units can provide a maximum of 500–750 ml/min of oxygen. But EverGo delivers 1,050 ml/min of oxygen through fixed bolus delivery, increasing the amount of oxygen with each setting. That means you have access to more oxygen with each breath – and a much lower chance of out-breathing the device.



\*Based on 20 breaths per minute and a typical device setting of 2.

