















# WHICH MASK IS THE BEST FOR YOU?

Mask Type	NASAL PILLOWS			NASAL MASKS						FULL FACE MASKS				
	Recommended when nasal bridge irritation occurs, or for claustrophobic users.			Suitable mask type for a majority of patients. Lightweight and comfortable, this mask type does not cover the whole face and promotes an effective breathing pattern with CPAP or Bi-Level						Recommended for patients with mouth leaks, blocked nasal passages and/or deviated septum				
Model	ResMed Swift FX	Fisher & Paykel Opus 360	Respironics OptiLife	ResMed Micro Mirage	ResMed Activa LT	ResMed Mirage SoftGel	Respironics ComfortGel	Respironics Easy Life	Fisher & Paykel Zest	ResMed Mirage Quattro	Fisher & Paykel Forma	Respironics ComfortGel Full Face	Respironics Full Life	Respironics Fit Life (Total Full Face)
														
Sizes	Xsmall, Small, Medium, Large	Small; Medium; Large	Petite; Small; Medium; Large	Small; Med/Large; Large-wide/XL	Small; Medium; Large-wide	Small; Medium; Large; Large Wide	Petite; Small; Medium; Large	Small; Medium; Medium/Wide; Large	Petite; Medium; Large	X-Small; Small; Medium; Large	Small; Medium; Large	Small; Medium; Large	Small; Medium; Large	Small; Large
Type of Cushions	Dual-wall nasal pillows	Built-in Soft-Lip pillows	Soft removable pillows and cradle cushion	Dual-wall cushion with Spring Air technology	ActiveCell Dual-wall cushion	DoubleGel Cushion	Customisable ComfortGel cushion	Dual-wall cushion	FlexiFit soft foam cushion	Dual-wall cushion with Spring Air technology	FlexiFit soft foam cushion	ComfortGel cushion	Dual textured cushion	Soft cushion
Elbow/ Tube/ Swivel	Interchangeable tube position	360° ball and socket elbow, interchangeable tube position	Interchangeable tube position	360° elbow with front attachment	360° elbow with front attachment	360° elbow with front attachment	360° swivel	Angled exhalation micro port make operation quieter and redirect air away from bed partner	360° elbow with front attachment	Simplified elbow with anti-asphyxia valve	360° elbow with front attachment	360° swivel	360° swivel allow to move freely	Build in radial diffusion leak ports deflect exhaled air away from bed partner
Forehead Support	Open face	Open face	Open face	MicroFit dial with 24 positions	MicroFit dial with 24 positions	MicroFit dial with 24 positions	Adjustable forehead support arm	Forehead pads adjust automatically, required minimum force to fit	Fit & Forget T-piece with 2 forehead cushions	MicroFit dial with 24 positions	Fit & Forget T-piece with 2 forehead cushions	Adjustable forehead support arm	None	None (Total full face mask)
Headgear	Contouring silicon headgear	Three-dimensional support structures	Two-tone fabric headgear, two strap slots & chin support band to support seal	Breath-O-Prene headgear, 4 adjustment points with 2 quick release clips	Breath-O-Prene headgear, 4 adjustment points with 2 quick release clips	Breath-O-Prene headgear, 4 adjustment points with 2 quick release clips	Soft flexible headgear straps improve comfort for all sleepers	Stretchgear headgear with Velcro straps	Stretchgear headgear with glider mechanism & Velcro straps	Set & forget headgear, 4 Velcro fasteners with 2 clips	Stretchgear headgear with glider mechanism & Velcro straps	Soft flexible headgear straps improve comfort for all sleepers	Soft Flexible headgear improve comfort for all sleepers	Soft Flexible headgear improve comfort for all sleepers

Washing Care

ALL: Hand wash in warm water with mild soap or dishwashing detergent. Rinse thoroughly and air dry.

## Call 1300 36 02 02



for the nearest AIR LIQUIDE *Healthcare* location or visit: [www.airliquidehealthcare.com.au](http://www.airliquidehealthcare.com.au)