







WHICH CPAP IS THE BEST FOR YOU?

Unit		Therapy Needs		Comfort Features		Convenience		Physical Characteristics		Humidifier	
Model		Type of PAP	Clinical Data Collection	Exhalations Comfort	Leak Compensation/Alert	Remote Update of Pressure Setting	Patient Reminder System	Size and Weight	Sound Level	Humidifier	Humidifier Chamber Features
Fisher & Paykel ICON Auto \$2,295		Auto Adjusting with Integrated Humidifier	F&P InfoSmart software via SmartStick	SensAwake pressure relief technology	Leak compensation	Yes with SmartStick	None	(H)x(W)x(D) 160mmX170mmX220mm 2.2kg with humidifier(420ml)	27dB	Integrated	ThermoSmart Technology included with heated breathing tube provides optimal humidity and reduced condensation
Fisher & Paykel ICON Premo \$1,450		Fixed pressure with Integrated Humidifier	F&P InfoSmart software via SmartStick	None	Leak compensation	Yes with SmartStick	None	(H)x(W)x(D) 160mmX170mmX220mm 2.2kg with humidifier(420ml)	27dB	Integrated	ThermoSmart Technology included with heated breathing tube provides optimal humidity and reduced condensation
Philips Respironics REMstar A-Flex \$1,995		Auto-adjusting	Encore Pro 2.1 via SD card reader	A-Flex	Leak compensation + mask off alert	Yes, with SD data card	Mask & filter replacement + service due date + supplier call	(H)x(W)x(D) 100mm x 140mm x 180mm; 0.9kg (does not include humidifier – 325mL)	25dB	System 1 \$300	Humidity control eliminating rainout; automatically adjusts to room temperature; complete water ingress protection
Philips Respironics REMstar Pro \$1,100		Fixed Pressure	Encore Pro 2.1 via SD card reader	C-Flex	Leak compensation + mask off alert	Yes, with SD data card	Mask & filter replacement + service due date + supplier call	(H)x(W)x(D) 100mm x 140mm x 180mm; 0.9kg (does not include humidifier – 325mL)	25dB	System 1 \$300	Humidity control eliminating rainout; automatically adjusts to room temperature; complete water ingress protection
ResMed S9 AutoSet \$2,099		Auto-adjusting	ResScan 3.10 software via SD card	EPR with Easy-Breathe technology	Leak compensation + mask off alert	Yes with SD data card	Mask & filter replacement + service due date + supplier call	(H)x(W)x(D) 86 mm x 140mm x 153mm; 0.81kg (does not include humidifier – 380mL)	25dB	H5i \$349	Climate control technology adjusts humidification to room temperature. Heated ClimateLine tube eliminate rainout
ResMed S9 Elite \$1,199		Fixed Pressure	ResScan 3.10 software via SD card	EPR with Easy-Breathe technology	Leak compensation + mask off alert	Yes with SD data card	Mask & filter replacement + service due date + supplier call	(H)x(W)x(D) 86 mm x 140mm x 153mm; 0.81kg (does not include humidifier – 380mL)	25dB	H5i \$349	Climate control technology adjusts humidification to room temperature. Heated ClimateLine tube eliminate rainout

▶ Start your Therapy Just **\$250** ◀



Call 1300 36 02 02

for the nearest AIR LIQUIDE Healthcare location or visit: www.airliquidehealthcare.com.au

