



Resuscitation

OxyAL™ Lifesaver



User-friendly resuscitation kit for swift response in an emergency

OxyAL™ Lifesaver is a compact resuscitation kit especially designed by Air Liquide Healthcare for quick emergency response.

- ✓ Compact and easy to carry
- ✓ Handy storage for accessories
- ✓ Available with reusable or disposable resuscitator kit



OxyAL™ with reusable resuscitator kit (Product code:3463SK) and OxyAL™ with disposable resuscitator kit (Product code: 3468SK), includes the following:

Item	Description	Product Code	Reusable/ Disposable
Resuscitator	<ul style="list-style-type: none"> • Reusable adult resuscitator with 1500ml ventilation bag, 60cmH₂O pop-off valve, size 5 silicone reusable mask, 2500ml reservoir bag & 2.1m non-kink tubing 	2857S	Reusable
	<ul style="list-style-type: none"> • Disposable adult resuscitator with 1500ml ventilation bag, 60cmH₂O pop-off valve, size 5 disposable mask, 2500ml reservoir bag & 2.1m non-kink tubing 	2858S	Disposable
OxyAL™ carry bag	<ul style="list-style-type: none"> • Compact, soft carry pack with shoulder straps and carry handles • Waterproof and easily cleaned 	1903S	Reusable
Guedel airways; size 0, 2, 4	<ul style="list-style-type: none"> • Airway adjunct used to maintain or open a patient's airway 	1893S 1894S 1895S	Disposable
Medium concentration mask	<ul style="list-style-type: none"> • Built for comfort and effective oxygen delivery • Threaded connectors eliminate the need to replace christmas trees on flowmeter 	154S	Disposable
Oxygen tube, 7 ft	<ul style="list-style-type: none"> • Three-channel safety tubing 		Disposable

Add-ons and accessories



Suction unit

Duo-O-Vac Suction Unit

Product code: 1205S

Duo-O-Vac Disposable Canister

Product code: 5239S

Twin-O-Vac Suction Unit

Product code: 1204S

Twin-O-Vac Reusable Canister

Product code: 2117S

Medium concentration mask

Disposable:

Size - paediatric

Product code: 157S

Size - adult

Product code: 154S



LEVEL 4, SUITES 4, 247 COWARD STREET
MASCOT NSW 2020

PH: 1300 36 02 02

W: airliquidehealthcare.com.au

E: alhenquiries@airliquide.com